



FACA-401-ISE

Fully Automatic Chemistry Analyzer



Labomed introduces an automatic biochemistry analyzer with varying testing capacities. FACA-401-ISE conducts 600 tests/hr (400 tests/hr + 200 ISE tests/hr, with minimum three ISE test (Na, K, Cl)) with True Random Access .


Features

- *Random Access and Direct Reading System.*
- *Automatic Washing System for Both Sample and Reagent Probe.*
- *Liquid Level Detection and Collision Protection for Both Sample and Reagent Probe.*
- *Reliable Cooling System for Reagent Chamber.*
- *Works with Windows XP and Windows 7*
- *Powerful Reaction Cuvette Washing System.*
- *User-Friendly Interface.*
- *ISE module and barcode scanner (optional)*
- *Enzyme Kinetic GPT and GOT testing*

Labomed, Inc. is certified by ISO 9001-2013, has CE Conformity and is FDA Licensed

Accessories

- *15 Reagent bottle 2 with holder*
- *80 Reaction cuvettes*
- *1000 Serum cups*
- *1 Tube*
- *4/4 Fuse (4A/8A)*
- *1 Lamp (13.8V/50W)*
- *2 Power Cables*
- *1 Com Cable*
- *1 Piston (1 ml)*
- *1 Bottle*
- *1 User Manual*
- *1 Operation Software CD*
- *1 Tool*
- *1 Sensor*
- *2 Detect Probes*
- *1 Sample Probe*
- *1 Black Cup*
- *1 Sinkage Block*



1 Reagent Tray

40 reagent positions for FACA-301 and FACA-261
 80 reagent positions for FACA-401
 Cooling system may work independently for 24 hours.

2 Sample Tray

60 sample positions for FACA-301 and FACA-261
 88 sample positions for FACA-401

3 Cuvette

Either high quality UV-transmitted plastic
 cuvettes or quartz glass cuvettes available.

4 Independent Mixer

1 mixer for FACA-301 and FACA-261
 2 mixers for FACA-401

5 Sample probe/ Reagent probe


1 sample probe and 1 reagent probe for FACA-301
 1 sample probe and 2 reagent probes for FACA-401
 Only 1 sample / reagent probe for FACA-261

6 ELE Module

ELE module available on request.

Cuvette Washing Station

8 channels automatic washing station.



Independent Mixer

Independent mixer ensuring quick and complete reaction.



Reaction Cuvette

High quality UV-transmitted plastic cuvettes, or quartz glass cuvette (optional).



Sample Probe / Reagent Probe

With liquid level sensor, anti-collision sensor and automatic cleaning function.




Reagent Bottle

More convenient and economical ready to use reagent kits.



Reagent Tray


With cooling reagent chamber.




Software




User-Friendly Interface
Compatible with Windows XP, Vista or Windows 7



Variable reading function




Carry Over Control Function




Automatic Alarm for Malfunctions

Automatic Washing Station




An 8-channel washing station ensures complete washing for cuvettes.

Environmental Protection Design



The waste is automatically classified and collected into different containers.

Automatic cuvette washing and blanking makes more precise results




Wavelength	Reading	Wash1	Wash2	Wash3	Wash4	Wash5	Wash6	Wash7	Wash8	Dry	Dry	New reading
340nm	1.3100	0.025	0.0092	0.0037	0.0022	0.0008	0.0006	0.0005	0.0002	0.0800	0.0850	/
405nm	0.6920	0.0112	0.0075	0.0033	0.0019	0.0010	0.0005	0.0003	0.0002	0.0800	0.0805	/
450nm	0.5281	0.0108	0.0066	0.0031	0.0015	0.0009	0.0003	0.0003	0.0002	0.0806	0.0805	/
510nm	0.7005	0.0113	0.0055	0.0027	0.0012	0.0006	0.0002	0.0002	0.0001	0.0810	0.0811	/
546nm	0.5368	0.0184	0.0072	0.0043	0.0013	0.0008	0.0002	0.0002	0.0002	0.0840	0.0843	0.5090
578nm	0.0136	0.0135	0.0064	0.0031	0.0012	0.0006	0.0002	0.0002	0.0002	0.0800	0.0800	/
620nm	0.0189	0.0189	0.0054	0.0026	0.0010	0.0005	0.0002	0.0002	0.0002	0.0800	0.0800	/


Test O.D. value-Real reagent blank=True O.D.value

General Features:

- New concept and more friendly software.
- New techniques and spare parts are used for much better quality and function.
- Automatic washing station for both sample and reagent probe.
- Liquid level detection and collision protection for both sample and reagent probe.
- Reliable cooling system for reagent chamber.
- Powerful reaction cuvette washing system.
- ISE module and barcode scanner(optional).
- Lower waste and more environment friendly.
- New concept test control system available (on request).



Reaction Cure Visible on line.




Cuvette Washing Station:
8 channels automatic washing station.



Independent mixer ensuring full reaction of sample / reagent solution.



High quality UV-transmitted plastic cuvettes, or quartz glass cuvette (optional).



Quickly heating up and thermostable technique to ensure reaction solution reach $37\pm 0.5^{\circ}\text{C}$ within 160s and stable at $37\pm 0.2^{\circ}\text{C}$ more than 1200s.



Sample probe / reagent probe with liquid level sensor, anti-collision sensor and automatic washing function.



Reagent tray can be easily placed in cooling chamber.

Technical Specifications

Model	FACA401-ISE
Assay Model	Random access with direct reading system
Method	End point, Kinetic, Fixed time, Immunoturbidimetric, 1-2 reagent, Multistandard, reagent/serum blank, etc
Assay Item	80 assay ISE mould (Na ⁺ , K ⁺ , Cl ⁻ / Li ⁺) available on request
Throughput	600 test/h (400 tests/hr + 200 ISE tests/hr, min. three ISE test (Na, K, Cl))
Sample Position	80 Including standard, control and stat; either serum cup or primary tube available
Sample Volume	1~100ul, 0.1 ul / step
Reagent Position	2 x 40
Reagent Volume	R1: 1~400ul, 1ul / step, R2: 1~400ul, 1ul / step
Reagent / Sample Probe	With liquid level sensor, collision sensor, teflon coating and automatic washing station
Reaction Time	0~999 seconds
Cuvette	High quality uv-transmitted plastic cuvette. Quartz glass cuvettes on request
Cuvette Washing	8 channels automatic washing system for reaction cuvettes
Probe And Mixer Washing	Independent washing station for each sample/reagent probe and mixer
Alarm	Automatic alarm for malfunction (eg: collision, reagent / sample insufficiency, full waste container, etc)
Calibration	Linear, Non-linear, Multi standard, K factor
Sample Dilution/Retest	Sample can be diluted and retested if results are out of range or sample is insufficient. Also with calibrator auto diluting function
Wavelength	340~810nm
Light Source	Halogen lamp
Absorbance Resolution	0.0001A
Absorbance Accuracy	±0.0003A(0~2.5A)
Repeatability (CV %)	≤2%
Quality Control	Levy Jennings quality control program at 3 levels
Bar Code	Optional
Printer	Multiple report formats available
Test Control Function	New card control system available on request
Language	English, French, Spanish, German, Russian, etc
System	Windows XP or Vista ; Windows 7
Power Supply	220V / 50Hz or 110V / 60Hz (optional)
Dimension	100×85×70 (cm)
Weight	137Kg