



LABOMED, INC.

www.labomed.com
spectro@labomed.com

LB-110 Monocular Biological Microscope

Introduction

With sharp image, competitive and reasonable unit price, LB-110 monocular biological microscopes are ideal instruments for student use. These microscopes are mainly designed for elementary school.

Specifications

Eyepiece:	Huygens Eyepiece H5x, H10x, H12.5x
Objective:	Achromatic Objective 10x, 40x
Nosepiece:	Triple Nosepiece
Stage:	Plain Stage with Slide Clips 110-120mm
Condenser:	Single Lens with Disc Diaphragm
Illumination System:	Plan-concave Mirror

