



www.labomed.com
spectro@labomed.com

Introduction

LB-116 Biological Monocular Microscopes are more comfortable and reliable operation with modern design. These microscopes are economical, with all basic features for educational use in elementary and high school.

Applications

These microscopes are economical microscope with all basic features for educational use in elementary and high school.

Technical Specifications

Viewing Head:	Monocular, Inclined at 45°, 360° Rotatable
Nosepiece:	Triple Nosepiece
Eyeiece:	Wide Field Eyepiece WF10×/18
Achromatic Objective:	Achromatic Objective 4×, 10×, 40×(S)
Stage:	Double Layers Mechanical Stage 110×120mm/ 70×20mm
Focusing:	Coaxial Coarse and Fine Adjustment
Condenser:	Abbe NA 1.25 with Iris Diaphragm & Blue Filter, Rack and Pinion Adjustment for Condenser
Illumination:	LED Illumination, Brightness Adjustable
Spare Parts:	Dust Cover
Package:	Styrofoam & Carton, Dimension 28×20×40cm, 3kg

