




# LABOMED, INC.

www.labomed.com  
spectro@labomed.com



## LB-210 Biological Binocular Educational Microscope with LED Illumination and Finite Optical System (Anti-Fungus)

Sliding Binocular Head, Inclined at 45°, 360° Rotatable

Built-in Revolving Quadruple Nosepiece.

Built-in Rotating Head

Optical Coating

Brass Gears

Wide Field Eyepiece WF10×/18, 4 Hole Nosepiece

Achromatic Objectives 4×, 10×, 40×(S), 100×

Double Layers Mechanical Stage 110x120mm/70X20mm Coaxial Coarse and Fine Adjustment

Abbe NA 1.25 with Iris Diaphragm & Blue Filter, Rack and Pinion Adjustment for Condenser

LED Illumination, Brightness Adjustable for Long Life and Energy Saving

Dust Cover

Styrofoam & Carton, Dimension 28x20x40cm, 3kg