




LABOMED, INC.

www.labomed.com
spectro@labomed.com



LB-275 Trinocular LED Fluorescent Biological Microscope with Wide Field and Infinite Plan Optical System (Infinity Color Corrected System) (Anti-Fungus)

Excellent fluorescent image with infinitive optical system.

LED illumination with 50,000 hours life span.

No need for cooling system and protection device due to LED cold light source.

465-476nm excited wavelength with high SNR.

Convenient to use, no need warm-up.

More fluorescent filters are available for different applications.