



LB-327

Zoom Stereo Microscope with Trinocular Head, WF10x Eyepiece, Zoom Objective and Universal Stand



LB-327 Zoom Stereo Microscope with Trinocular Head, WF10x Eyepiece, Zoom Objective and Universal Stand

LB-327 Zoom Stereo Microscope with Trinocular Head, WF10x Eyepiece, Zoom Objective and Universal Stand offer sharp 3D images which are very clear throughout the zoom range. These microscopes are very popular and cost effective. Optional eyepieces and auxiliary objectives can expand magnification range and working distances. Cold light and ring light can be chosen for this microscope.

APPLICATION

LB-327 Zoom Stereo Microscope with Trinocular Head, WF10x Eyepiece, Zoom Objective and Universal Stand are widely used in education, lab research, biology, metallurgy, engineering, chemistry, manufacturing, and in the medical, forensic science and veterinary industries. The microscopes also can be used for circuit board repair and inspection, SMT work, electronics inspection, dissection, coin collecting, gemology and gemstone setting, engraving, mobile phone repair and inspection of small parts.

FEATURES

1. 7×-45× zoom magnification power with sharp images, can be extended to 3.5×-270× with optional eyepiece and auxiliary objective.
2. High eyepoint WF10×/20mm eyepiece.
3. 100mm Long working distance.
4. Ergonomic design, sharp image, wide viewing field, high depth of field and easy to operate.
5. Ideal instrument in education, medical and industrial field.

SPECIFICATION

Item	Specification
Viewing Head	Trinocular head, inclined at 45°, Interpupillary Distance, 54-76mm, 2:8, diopter adjustment for both tubes, 30mm tube
Eyepiece	WF10×/ 20mm eyepiece (micrometer is optional)
Objective	Zoom objective 0.7×-4.5×
Zoom Ratio	01:06.3
Working Distance	100mm
Head Mount	76mm
Focusing Arm	Coarse focusing, focusing range 50mm
Pillar Stand	ST1 single arm universal stand, no illumination, can expand the operating range and observe objects in any direction.
	Vertical column height: 380mm
	Column diameter: φ32mm
	Horizontal movement: 235mm
	Base: 230x230x40mm
	With FA1 coarse focusing arm(76mm mounting size)

Package	1set/1carton, 49.5cm*48.5cm*29cm, Net/Gross Weight: 19/20kg
---------	--

SAMPLE IMAGES

