




LABOMED, INC.

www.labomed.com
spectro@labomed.com



LB-341 Trinocular Zoom Stereo Microscope

With high quality optical components, providing sharp and striking contrast image and ensure flat image at large depth of field. Ease the tiredness with ergonomic design.

With LED light for both incident and transmitted illumination, providing even illumination and life expectancy of the LED light can reach 6000 hours.

With Complete eyepieces, objectives and accessories, it is an ideal instrument for industry assembling, inspection and teaching field.

For accessories, please go to LB-150 Cold Light Illumination