




www.labomed.com
spectro@labomed.com

5MP Digital Camera for LB-356 Digital Zoom Stereo Microscope

1. Real-Time Measurement Software

Equipped with international standard C-type interface, it can be conveniently thread on both C-mount coupler and eyepiece adapter. It is compatible with the monocular, binocular, trinocular and any kind of microscopes, such as stereo, video, dissecting microscope or telescopes.

- ① Preview Image – show live video of the camera.
- ② Tool Panel – Control panels, capture, adjust parameter, measurement, etc.
- ③ Thumbnail Panel – Show captured images and videos, touch it to browse the files.
- ④ System Tool – Switch Camera, full screen, configure.



2. Particle Analysis Software

Touch the "YWImage" icon to open the particle analysis application.

