




LABOMED, INC.

www.labomed.com  
spectro@labomed.com



## LB-372 Wide Field Digital Zoom Stereo Microscope with Universal Stand

- With high quality optical components, providing sharp and high contrast image at super-long working distance, super-large viewing field, large depth of focus, reduce the fatigue and improve the working efficiency.
- With LED light, Providing even illumination and life expectancy can reach 60000 hours.
- Various objectives are available, 2×, 4×, 6× objectives are standard, 8×, 10×, 15× and 6×SL objectives are optional. The objectives can be conveniently changed.