




LABOMED, INC.

www.labomed.com
spectro@labomed.com



FOR NON-TRANSPARENT SPECIMENS

LB-592 Polarizing Trinocular Microscope with Infinite Optical System

Trinocular head, digital cameras can be easily used on this microscope for image and video taking.
Supply high quality and contrast image with special strain free objectives.

Can be used for single polarization, orthogonal polarization and conoscopic observation.

Transmitted illumination, 24V/100W long life halogen lamp.

Stage with 45° located, attachable mechanical stage with large moving range, no interference when change objectives.

Gypsum λ plate, mica $\lambda/4$ plate, quartz wedge plate are equipped.