




# LABOMED, INC.

www.labomed.com  
spectro@labomed.com



## LB-612 Binocular Metallurgical Microscope with Infinite Plan Optical System, Extra Wide Field, 5 objectives, Backward Quintuple Nosepiece, Transmitted and Reflected Kohler and Halogen Illumination, Polarizer and Analyzer

Infinite optical system provides excellent optical functions.

Laboratory metallurgical microscope, including bright field, dark field, polarization observation system.

Powerful transmitted and reflected system with Kohler illumination.

Ideal instrument for industry inspection and science research.