

# Labomed, Inc. 2019 Price List

## Cuvette Price List



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price (USD)		Retail Price (USD)	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G1	Q1	45x12.5x3.5	1	10	0.35	23	79	46	158
G2	Q2	45x12.5x4.5	2	10	0.70	23	79	46	158
G3	Q3	45x12.5x7.5	5	10	1.70	23	79	46	158
G4	Q4	45x12.5x12.5	10	10	3.50	17	45	34	90
G5	Q5	45x12.5x22.5	20	10	7.00	29	69	58	138
G6	Q6	45x12.5x32.5	30	10	10.50	29	81	58	162
G7	Q7	45x12.5x42.5	40	10	14.00	35	92	70	184
G8	Q8	45x12.5x52.5	50	10	17.50	52	110	104	220
G9	Q9	45x12.5x102.5	100	10	35.00	63	224	126	448

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G13	Q13	45x12.5x7.5	5	10	1.70	46	102	92	204
G14	Q14	45x12.5x12.5	10	10	3.50	46	102	92	204
G15	Q15	45x12.5x22.5	20	10	7.00	53	114	106	228
G16	Q16	45x12.5x32.5	30	10	10.50	58	127	116	254
G17	Q17	45x12.5x42.5	40	10	14.00	60	138	120	276
G18	Q18	45x12.5x52.5	50	10	17.50	67	161	134	322
G19	Q19	45x12.5x102.5	100	10	35.00	138	380	276	760

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm




Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G23	Q23	45x12.5x7.5	5	4	0.70	52	63	104	126
G24	Q24	45x12.5x12.5	10	4	1.40	52	63	104	126
G25	Q25	45x12.5x22.5	20	4	2.80	58	86	116	172
G26	Q26	45x12.5x32.5	30	4	4.20	69	98	138	196
G27	Q27	45x12.5x42.5	40	4	5.60	63	99	126	198
G28	Q28	45x12.5x52.5	50	4	7.00	69	114	138	228
G29	Q29	45x12.5x102.5	100	4	14.00	150	253	300	506

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G33	Q33	45x12.5x7.5	5	4	0.70	52	102	104	204
G34	Q34	45x12.5x12.5	10	4	1.40	52	102	104	204
G35	Q35	45x12.5x22.5	20	4	2.80	58	114	116	228
G36	Q36	45x12.5x32.5	30	4	4.20	63	127	126	254
G37	Q37	45x12.5x42.5	40	4	5.60	69	138	138	276
G38	Q38	45x12.5x52.5	50	4	7.00	79	161	158	322
G39	Q39	45x12.5x102.5	100	4	14.00	138	380	276	760

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G43	Q43	45x12.5x7.5	5	2	0.35	46	69	92	138
G44	Q44	45x12.5x12.5	10	2	0.70	54	69	108	138
G45	Q45	45x12.5x22.5	20	2	1.40	58	81	116	162
G46	Q46	45x12.5x32.5	30	2	2.10	58	92	116	184
G47	Q47	45x12.5x42.5	40	2	2.80	58	99	116	198
G48	Q48	45x12.5x52.5	50	2	3.50	69	127	138	254

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm

Micro cell with black walls &amp; teflon stopper



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G53	Q53	45x12.5x7.5	5	2	0.35	58	102	116	204
G54	Q54	45x12.5x12.5	10	2	0.70	69	102	138	204
G55	Q55	45x12.5x22.5	20	2	1.40	75	114	150	228
G56	Q56	45x12.5x32.5	30	2	2.10	75	127	150	254
G57	Q57	45x12.5x42.5	40	2	2.80	75	138	150	276
G58	Q58	45x12.5x52.5	50	2	3.50	75	161	150	322
G59	Q59	45x12.5x102.5	100	2	7.00	127	380	254	760

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm

Self masking continuous flowthrough cell



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G63	Q63	45x12.5x7.5	5	3	0.035	69	81	138	162
G64	Q64	45x12.5x12.5	10	3	0.07	69	92	138	184
G65	Q65	45x12.5x22.5	20	3	0.14	69	78	138	156
G66	Q66	45x12.5x32.5	30	3	0.21	75	83	150	166
G67	Q67	45x12.5x7.5	5	2	0.015	81	81	162	162
G68	Q68	45x12.5x12.5	10	2	0.03	69	75	138	150
G69	Q69	45x12.5x22.5	20	2	0.06	81	92	162	184
G70	Q70	45x12.5x12.5	10	1.5	0.018	81	92	162	184

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm

Self masking continuous flowthrough cell



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G73	Q73	45x12.5x7.5	5	4x12	0.24	69	81	138	162
G74	Q74	45x12.5x12.5	10	4x12	0.48	69	92	138	184
G75	Q75	45x12.5x22.5	20	4x12	1.00	69	92	138	184
G76	Q76	45x12.5x32.5	30	4x12	1.50	81	91	162	182
G79	Q79	45x12.5x7.5	5	2x12	0.12	69	81	138	162
G80	Q80	45x12.5x12.5	10	2x12	0.24	69	81	138	162

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm

Cylindrical cells



Item No.		Outside Dimensions LxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G85	Q85	22.5x12.5	10	19	2.80	58	92	116	184
G86	Q86	22.5x22.5	20	19	5.70	69	113	138	226
G87	Q87	22.5x52.5	50	19	14.00	81	99	161	198
G88	Q88	22.5x102.5	100	19	28.00	121	138	242	276

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm

Cylindrical cells one stopper



Item No.		Outside Dimensions LxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G95	Q95	22.5x12.5	10	19	2.80	58	69	116	138
G96	Q96	22.5x22.5	20	19	5.70	58	81	116	162
G97	Q97	22.5x52.5	50	19	14.00	69	92	138	184
G98	Q98	22.5x102.5	100	19	28.00	121	150	242	300

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm

Standard cells with lid and round bottom



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G101	Q101	45x12.5x3.5	1	10	0.35	23	104	46	208
G102	Q102	45x12.5x4.5	2	10	0.70	23	75	46	150
G103	Q103	45x12.5x7.5	5	10	1.70	23	63	46	126
G104	Q104	45x12.5x12.5	10	10	3.50	23	45	46	90
G105	Q105	45x12.5x22.5	20	10	7.00	29	69	58	138
G106	Q106	45x12.5x32.5	30	10	10.50	29	81	58	162
G107	Q107	45x12.5x42.5	40	10	14.00	35	92	70	184

G108	Q108	45x12.5x52.5	50	10	17.50	52	110	104	220
G109	Q109	45x12.5x102.5	100	10	35.50	63	224	126	448

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm

Standard cell with Teflon stopper and round bottom



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G113	Q113	45x12.5x7.5	5	10	1.70	46	127	92	254
G114	Q114	45x12.5x12.5	10	10	3.50	46	127	92	254
G115	Q115	45x12.5x22.5	20	10	6.00	53	150	106	300
G116	Q116	45x12.5x32.5	30	10	10.50	58	166	116	332
G117	Q117	45x12.5x42.5	40	10	14.00	60	198	120	396
G118	Q118	45x12.5x52.5	50	10	17.50	67	225	134	450
G119	Q119	45x12.5x102.5	100	10	35.00	138	380	276	760

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm

Micro cell with black walls and with lid



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G123	Q123	45x12.5x7.5	5	1	0.17	46	81	92	162
G124	Q124	45x12.5x12.5	10	1	0.35	46	81	92	162
G125	Q125	45x12.5x22.5	20	1	0.70	58	92	116	184
G126	Q126	45x12.5x32.5	30	1	1.00	63	98	126	196
G127	Q127	45x12.5x42.5	40	1	1.40	69	110	138	220
G128	Q128	45x12.5x52.5	50	1	1.75	81	127	162	254
G129	Q129	45x12.5x102.5	100	1	3.50	150	276	300	552

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm

Micro cell with black walls and lid



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G133	Q133	45x12.5x7.5	5	3	0.50	46	69	92	138
G134	Q134	45x12.5x12.5	10	3	1.00	46	75	92	150
G135	Q135	45x12.5x22.5	20	3	2.10	58	81	116	162
G136	Q136	45x12.5x32.5	30	3	3.15	69	92	138	184
G137	Q137	45x12.5x42.5	40	3	4.20	75	104	150	208
G138	Q138	45x12.5x52.5	50	3	5.25	81	127	162	254
G139	Q139	45x12.5x102.5	100	3	10.50	127	276	254	552

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm

Semi I - Micro cell with black walls & lid



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G143	Q143	45x12.5x7.5	5	5	0.85	46	69	92	138
G144	Q144	45x12.5x12.5	10	5	1.75	46	69	92	138
G145	Q145	45x12.5x22.5	20	5	3.50	52	81	104	162
G146	Q146	45x12.5x32.5	30	5	5.25	58	92	116	184
G147	Q147	45x12.5x42.5	40	5	7.00	69	104	138	208
G148	Q148	45x12.5x52.5	50	5	8.75	81	115	162	230
G149	Q149	45x12.5x102.5	100	5	17.50	138	276	276	552

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm

Self masking continuous flowthrough cell



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G154	Q154	45x12.5x12.5	10	1	0.01	69	81	138	162
G155	Q155	45x12.5x22.5	20	1	0.016	69	92	138	184
G156	Q156	45x12.5x32.5	30	1	0.025	76	98	152	196
G164	Q164	45x12.5x12.5	10	1.5	0.018	63	92	126	184
G165	Q165	45x12.5x22.5	20	1.5	0.036	69	92	138	184
G166	Q166	45x12.5x32.5	30	1.5	0.055	81	92	162	184
G174	Q174	45x12.5x12.5	10	2	0.030	86	104	172	208
G175	Q175	45x12.5x22.5	20	2	0.060	92	109	184	218

G176	Q176	45x12.5x32.5	30	2	0.090	98	115	196	230
------	------	--------------	----	---	-------	----	-----	-----	-----

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm

Standard fluorimeter cell with lid



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G203	Q203	45x12.5x7.5	5	10	1.70	23	92	46	184
G204	Q204	45x12.5x12.5	10	10	3.50	35	98	70	196
G205	Q205	45x12.5x22.5	20	10	7.00	35	104	70	208
G206	Q206	45x12.5x32.5	30	10	10.50	46	127	92	254
G207	Q207	45x12.5x42.5	40	10	14.00	58	150	116	300
G208	Q208	45x12.5x52.5	50	10	17.50	69	173	138	346
G209	Q209	45x12.5x102.5	100	10	35.00	115	288	230	576

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm

Standard fluorimeter cell w/ teflon stopper



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G213	Q213	45x12.5x7.5	5	10	1.70	58	92	116	184
G214	Q214	45x12.5x12.5	10	10	3.50	58	92	116	184
G215	Q215	45x12.5x22.5	20	10	7.00	58	104	116	208
G216	Q216	45x12.5x32.5	30	10	10.50	69	138	138	276
G217	Q217	45x12.5x42.5	40	10	14.00	75	161	150	322
G218	Q218	45x12.5x52.5	50	10	17.50	81	184	162	368
G219	Q219	45x12.5x102.5	100	10	35.00	138	334	276	668

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm

Micro cell with frosted walls and lid



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G233	Q233	45x12.5x7.5	5	1	0.17	43	69	86	138
G234	Q234	45x12.5x12.5	10	1	0.35	54	89	108	178
G234.5	Q234.5	45x12.5x22.5	10	1	0.50	58	86	116	172
G235	Q235	45x12.5x22.5	20	1	0.70	46	92	92	184
G236	Q236	45x12.5x32.5	30	1	1.00	58	104	116	208
G237	Q237	45x12.5x42.5	40	1	1.40	58	109	116	218
G238	Q238	45x12.5x52.5	50	1	1.75	64	113	128	226
G239	Q239	45x12.5x102.5	100	1	3.50	127	265	254	530

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm

Micro cell with frosted walls and lid



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G243	Q243	45x12.5x7.5	5	2	0.35	46	75	92	150
G244	Q244	45x12.5x12.5	10	2	0.70	46	75	92	150
G245	Q245	45x12.5x22.5	20	2	1.40	58	81	116	162
G246	Q246	45x12.5x32.5	30	2	2.10	58	92	116	184
G247	Q247	45x12.5x42.5	40	2	2.80	69	104	138	208
G248	Q248	45x12.5x52.5	50	2	3.50	69	127	138	254
G249	Q249	45x12.5x102.5	100	2	7.00	138	276	276	552

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm

Micro cell with frosted walls and lid



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G253	Q253	45x12.5x7.5	5	3	0.50	46	75	92	150
G254	Q254	45x12.5x12.5	10	3	1.00	46	75	92	150
G255	Q255	45x12.5x22.5	20	3	2.10	58	81	116	162
G256	Q256	45x12.5x32.5	30	3	3.50	58	92	116	184
G257	Q257	45x12.5x42.5	40	3	4.20	69	104	138	208
G258	Q258	45x12.5x52.5	50	3	5.25	69	127	138	254
G259	Q259	45x12.5x102.5	100	3	10.50	138	276	276	552

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm

Micro cell with lid



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G263	Q263	45x12.5x7.5	5	4	0.70	46	75	92	150
G264	Q264	45x12.5x12.5	10	4	1.40	46	75	92	150
G265	Q265	45x12.5x22.5	20	4	2.80	58	81	116	162
G266	Q266	45x12.5x32.5	30	4	4.20	58	92	116	184
G267	Q267	45x12.5x42.5	40	4	5.60	69	104	138	208
G268	Q268	45x12.5x52.5	50	4	7.00	69	127	138	254
G269	Q269	45x12.5x102.5	100	4	14.00	138	276	276	552

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm

Semi Micro cell with frosted walls &amp; lid



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G273	Q273	45x12.5x7.5	5	5	0.85	46	75	92	150
G274	Q274	45x12.5x12.5	10	5	1.70	46	75	92	150
G275	Q275	45x12.5x22.5	20	5	3.50	58	81	115	162
G276	Q276	45x12.5x32.5	30	5	5.20	58	92	115	184
G277	Q277	45x12.5x42.5	40	5	7.00	69	104	138	208
G278	Q278	45x12.5x52.5	50	5	8.75	69	127	138	254
G279	Q279	45x12.5x102.5	100	5	17.50	138	276	276	552

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm

Semi Micro cell with frosted walls &amp; teflon stopper



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G283	Q283	45x12.5x7.5	5	2	0.35	58	102	116	204
G284	Q284	45x12.5x12.5	10	2	0.70	58	102	116	204
G285	Q285	45x12.5x22.5	20	2	1.40	63	114	126	228
G286	Q286	45x12.5x32.5	30	2	2.10	69	127	138	254
G287	Q287	45x12.5x42.5	40	2	2.80	75	138	150	276
G288	Q288	45x12.5x52.5	50	2	3.50	81	161	162	322
G289	Q289	45x12.5x102.5	100	2	7.00	127	380	254	760

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm

Semi Micro cell with frosted walls &amp; teflon stopper



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G293	Q293	45x12.5x7.5	5	4	0.75	56	102	112	204
G294	Q294	45x12.5x12.5	10	4	1.40	56	102	112	204
G295	Q295	45x12.5x22.5	20	4	2.80	63	114	126	228
G296	Q296	45x12.5x32.5	30	4	4.20	69	127	138	254
G297	Q297	45x12.5x42.5	40	4	5.60	75	138	150	276
G298	Q298	45x12.5x52.5	50	4	7.00	81	161	162	322
G299	Q299	45x12.5x102.5	100	4	14.00	144	380	288	760

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm

Large cell



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G321	Q321	40x20x10.4	6.35	16	3.50	53	150	106	300
G322	Q322	40x20x16.7	12.7	16	7.10	53	150	106	300
G323	Q323	40x20x29.4	25.4	16	14.20	58	161	116	322
G324	Q324	40x20x54.8	50.8	16	28.40	58	161	116	322
G325	Q325	40x20x137.4	133.5	16	74.60	81	224	162	448
G326	Q326	40x20x14	10	16	5.60	58	224	116	448
G327	Q327	40x20x104	100	16	56.00	81	224	162	448

Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm

Micro cell with frosted walls &amp; teflon stopper



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G373	Q373	45x12.5x7.5	5	1	0.17	56	102	112	204
G374	Q374	45x12.5x12.5	10	1	0.35	56	102	112	204
G375	Q375	45x12.5x22.5	20	1	0.70	63	114	126	228
G376	Q376	45x12.5x32.5	30	1	1.00	69	127	138	254
G377	Q377	45x12.5x42.5	40	1	1.40	75	138	150	276
G378	Q378	45x12.5x52.5	50	1	1.75	81	161	162	322
G379	Q379	45x12.5x102.5	100	1	3.50	144	380	288	760

**Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm**

Micro cell with frosted walls &amp; teflon stopper



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G383	Q383	45x12.5x7.5	5	3	0.50	56	102	112	204
G384	Q384	45x12.5x12.5	10	3	1.00	56	102	112	204
G385	Q385	45x12.5x22.5	20	3	2.10	63	114	126	228
G386	Q386	45x12.5x32.5	30	3	3.15	69	127	138	254
G387	Q387	45x12.5x42.5	40	3	4.20	75	138	150	276
G388	Q388	45x12.5x52.5	50	3	5.25	81	161	162	322
G389	Q389	45x12.5x102.5	100	3	10.50	144	380	288	760

**Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm**

Micro cell with frosted walls &amp; teflon stopper



Item No.		Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price	
Glass	Quartz					Glass	Quartz	Glass	Quartz
G393	Q393	45x12.5x7.5	5	5	0.85	56	102	112	204
G394	Q394	45x12.5x12.5	10	5	1.70	56	102	112	204
G395	Q395	45x12.5x22.5	20	5	3.50	63	114	126	228
G396	Q396	45x12.5x32.5	30	5	5.20	69	127	138	254
G397	Q397	45x12.5x42.5	40	5	7.00	75	138	150	276
G398	Q398	45x12.5x52.5	50	5	8.75	81	161	162	322
G399	Q399	45x12.5x102.5	100	5	17.50	144	380	288	760

**Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm**

Micro cell with frosted walls and lid



Item No. G or Q	Outside Dimensions HxWxD mm	Path Length (mm)	Inside Width (mm)	Volume (mm)	Dist. Price		Retail Price		
					Glass	Quartz	Glass	Quartz	
Q-674-10	45x12.5x12.5	10	10	0.01	10.00	N/A	109	N/A	218
Q-674-20	45x12.5x12.5	10	10	0.02	20.00	N/A	109	N/A	218
Q-674-50	45x12.5x12.5	10	10	0.05	50.00	N/A	109	N/A	218
Q-674-100	45x12.5x12.5	10	10	0.10	100.00	N/A	109	N/A	218
Q-674-200	45x12.5x12.5	10	10	0.20	200.00	N/A	109	N/A	218
Q-674-500	45x12.5x12.5	10	10	0.50	500.00	56	92	112	184
Q-674-1000	45x12.5x12.5	10	10	1.00	1000.00	56	92	112	184

**Q = Quartz / Cells for 190-2500nm / G = Optical Glass / Cells for 320-2500nm**