



Features :

- Random access and direct reading system.
- Automatic washing system for both sample and reagent probe.
- Liquid level detection and collision protection for both sample and reagent probe.
- Reliable cooling system for reagent chamber.
- Powerful reaction cuvette washing system.
- User-friendly interface.
- ISE module and barcode scanner (optional).

• CE • ISO9001:2000

FACA-401

FACA-301

FACA-261

FULLY AUTOMATIC CHEMISTRY ANALYZER

Cuvette Washing Station

8 channels automatic washing station.



Independent Mixer

Independent mixer ensuring quick and complete reaction.




Reaction Cuvette

High quality UV-transmitted plastic cuvettes, or quartz glass cuvette (optional).



Sample Probe / Reagent Probe

With liquid level sensor, anti-collision sensor and automatic cleaning function.



Reagent Bottle

More convenient and economical ready to use reagent kits.



Reagent Tray

With cooling reagent chamber.

