


# LABOMED



## 1 Reagent Tray

40 reagent positions for DS-301 and DS-261.  
80 reagent positions for DS-401.  
Cooling system may work independently for 24 hours.

## 2 Sample Tray

60 sample positions for DS-301 and DS-261.  
88 sample positions for DS-401.

## 3 Cuvette

Either high quality UV-transmitted plastic  
cuvettes or quartz glass cuvettes available.

## 4 Independent Mixer

1 mixer for DS-301 and DS-261.  
2 mixers for DS-401.

## 5 Sample probe/ Reagent probe

1 sample probe and 1 reagent probe for DS-301.  
1 sample probe and 2 reagent probes for DS-401.  
Only 1 sample / reagent probe for DS-261.

## 6 ISE Module

ISE module available on request.