


H-2150 Fully Automatic 5-part Differential Hematology Analyzer



H-2150 Fully Automatic 5-part Differential Hematology Analyzer is a 5-part diff. multi-parameter automated Hematology Analyzer designed for clinical laboratories. The analyzer uses light scatter and impedance technology for WBC differential and impedance method for RBC and PLT. With the most updated technology and concepts, it best meets most clinical laboratory requirements for an automated hematology analyzer. The equipment offers user-friendly software, powerful functions, stable part differential, and accurate results.