




**LABOMED, INC.**

[www.labomed.com](http://www.labomed.com)  
[spectro@labomed.com](mailto:spectro@labomed.com)



## **LB-101 Monocular Zoom Microscope**

|                              |                                                                 |
|------------------------------|-----------------------------------------------------------------|
| <b>Eyepiece:</b>             | <b>Eyepiece WF10x</b>                                           |
| <b>Zoom Ratio:</b>           | <b>Zoom Ratio: 1:6.5 (0.7×-4.5×)</b>                            |
| <b>Total Magnification:</b>  | <b>7×-45×</b>                                                   |
| <b>Working Distance:</b>     | <b>90mm</b>                                                     |
| <b>Stand Focusing Range:</b> | <b>55mm</b>                                                     |
| <b>Stand Clips:</b>          | <b>90mm Frosted Glass or Black/ White Plate &amp; Clips</b>     |
| <b>Illumination:</b>         | <b>Reflect Illumination, 6V/15W Halogen Lamp with Reflector</b> |
| <b>Bottom Illumination:</b>  | <b>6V/20W Halogen Lamp</b>                                      |
| <b>Video System:</b>         | <b>C Mount</b>                                                  |