





LABOMED, INC.

www.labomed.com
spectro@labomed.com



LB-1020 Digital Comparison Microscope

LB-1020 Digital Comparison Microscope is a high level Digital Comparison Microscope, it has a built-in digital image analysis system, it is basically used for forensic science and police departments. It is also provided to these departments such as police training schools, banks, tax, print, coinage, archaeology and etc.

3 Megapixel CCD camera, supports USB/AV output, have camera connector

2.5×-130×, LCD display

24V 300W cold light source, adjustable

This microscope comes with camera and software and illuminates by R, G, B LED cold light, and can change color. This microscope focuses electronically. It focuses by an electrical button. The objectives are changed, one by one.