




LABOMED, INC.

www.labomed.com
spectro@labomed.com



LB-113 Monocular Biological Microscope

Viewing Head:	Monocular, Inclined at 45°, 360° Rotatable
Nosepiece:	Triple Nosepiece
Eyepiece:	Wide Field Eyepiece WF10×/18
Achromatic Objective:	Achromatic Objective 4×, 10×, 40×(S)
Stage:	Plain Stage with Slide Clips 110×120mm
Focusing:	Coaxial Coarse and Fine Adjustment, Moving Range 10mm
Condenser:	Single Lens NA 0.65 with Disc Diaphragm
Illumination:	LED Illumination, Brightness Adjustable
Spare Parts:	Dust Cover
Package:	Styrofoam & Carton, Dimension 28×20×40cm, 3kg