




LABOMED, INC.

www.labomed.com
spectro@labomed.com



Swivel Head LCD



LB-1250 LCD Digital Biological Microscope

Variety observation mode: Observe image on LCD screen and ocular at the same time; Support USB, AV, S-Video and Audio output and display real-time image on TV and PC at the same time.

2.5" TFT LCD screen with 4× digital zooming, rotatable at maximum 100° elevation angle and 180° pivoting.

Built-in 3M pixel digital camera facilitates capturing images, video and audio with no messy wires or power cords. Expandable to Max 1G SD Card slot for more storage capacity.

Support Windows 2000/XP/Vista/Windows 7 and 8 and equipped with powerful image processing software. Integrated traditional optics and digital technology.