




LABOMED, INC.

www.labomed.com
spectro@labomed.com



Swivel Head LCD



LB-1251 Compound LCD Digital Biological Microscope (3.0MP) with Infinite Semi-Plan Optical System (Infinity Color Corrected System) and Wide Field

Variety observation mode: Observe image on LCD screen and ocular at the same time; Support USB, AV, S-Video and Audio output and display real-time image on TV and PC at the same time.
2.5" TFT LCD screen with 4× digital zooming, rotatable at maximum 100° elevation angle and 180° pivoting.
Built-in 3M pixel digital camera facilitates capturing images, video and audio with no messy wires or power cords.
Expandable to Max 1G SD Card slot for more storage capacity.
Support Windows 2000/XP/Vista and equipped with powerful image processing software.
Integrated traditional optics and digital technology.