




# LABOMED, INC.

www.labomed.com  
spectro@labomed.com



## **LB-206 Binocular Digital Biological Microscope with Seidentopf Binocular Head, Wide Field Eyepiece, Backward Quadruple Nosepiece, Infinity E-Plan Objectives, and 3.2MP Camera**

Excellent infinite optical system, outstanding image quality with high resolution and definition.

Dustproof, moisture-proof, built-in chip.

Light distribution 20:80, eyepiece observation synchronized with digital imaging without switching.

Built-in 3.2MP digital camera system, 1/2" CMOS color sensor.

USB2.0 connection without external power supply.

Available for WINDOWS 2000/XP/7/Vista system.