




LABOMED, INC.

www.labomed.com
spectro@labomed.com



LB-237 Binocular Biological Digital Microscope with Wide Field, Semi-Plan Optical System and Camera and Software (1.3MP)

High resolution, real-time image display by 1.3 Megapixel CMOS digital camera connecting to PC through USB2.0.

Sharpe Image with Semi-plan Achromatic Objectives.

Easy to replace lamp.

High sensitive camera can be used for optional attachments such as dark field, phase contrast, Epi-fluorescence and so on.

Professional and powerful micro-image process and analysis software.