




LABOMED, INC.

www.labomed.com
spectro@labomed.com



LB-243 Trinocular Biological Microscope with Simple Phase Contrast Kit, Digital Camera (3.2MP) and Infinite Optical System

Extra wide Field Eyepiece EW10x/18 with Diopter Adjustment as Optional
Comfortable Observation with Eyepoint and interpupillary adjustable
Sliding-in Centerable Condenser Adjustable
Handle for easy carrying
Lamp replacement by directly removing housing
The Simple Phase Contrast Kit
The Color Digital Camera is 3.2MP.