




LABOMED, INC.

www.labomed.com
spectro@labomed.com

LB-260 Binocular Medical Microscope

Optional Features

	<p>OPTIONAL:</p> <p>Monocular, Binocular and Trinocular Head Available for Different Requirements.</p>
	<p>Big Operating Space With Backward Nosepiece. Easier Comparing Test with Large Size Double Layer Mechanical Stage.</p>
	<p>More Comfortable Operation with Low Position CA and FA Knob, Low Stage Control and Narrow Base. Even Illumination with Modern Design, Kohler Illumination and Aspheric Collector can be Adapted.</p>
	<p>Aluminum Carrying Case Available for High Level Packing Requirement</p>