



www.labomed.com
spectro@labomed.com

LB-931 Multi-Head Dual View (Front and Back) Microscope (Anti-Fungus)

Introduction

LB-931 Multi-Head Dual View (Front and Back) Microscopes (Anti-Fungus) are equipped with multi-heads for more persons observing at the same time. The optical system is high quality and reliable.

Applications

LB-931 Multi-Head Dual View (Front and Back) Microscopes (Anti-Fungus) are mainly used in medical teaching and biological teaching. They also can be used for biological analysis for more experts at the same time. They can be widely used in medical and sanitary establishments, laboratories, institutes, colleges and universities.

Technical Specifications

Viewing Head:	Sliding Binocular Head, inclined at 45°
Eyepiece:	Wide Field Eyepiece WF10×
Nosepiece:	Quadruple Nosepiece
Objective:	Achromatic Objective 4×, 10×, 40×, 100×
Condenser:	Abbe Condenser NA1.25
Focusing:	Coaxial Coarse & Fine Focus Adjustment, Fine Division 0.002mm
Stage:	Double Layers Mechanical Stage 140×140mm/ 75×50mm Plain Stage 140×140mm/ 75×50mm
Illumination:	LED Halogen Lamp 6V/ 20W (optional)

