

No.	Description	Abbreviation	Method	New package			Storage		Sale Price
				R1	R2	STD	°C = Temperature of Storage	Lifetime	
				x= Quantity in box R1 and R2 should be mixed before testing					
	Liver Function			R1	R2	REFERENCE	°C	Lifetime	
	Alanine Aminotransferase	ALT/GPT	IFCC Method	32ml x4	8ml x4		2-8	1year	\$18.00
	Aspartate Aminotransferase	AST/GOT	IFCC Method	32ml x4	8ml x4		2-8	1year	\$18.00
	γ-Glutamyltransferase	γ-GT/GGT	IFCC Method	32ml x4	8ml x4		2-8	1year	\$40.00
	Alkaline Phosphatase	ALP/AKP	IFCC Method	32ml x4	8ml x4		2-8	1year	\$24.00
	Total Bilirubin*	TBILI	Diazo Salt Method	35ml x4	7ml x4	2ml x1	2-8	2year	\$20.00
	Total Bilirubin*	TBILI	Vanadate Method	32ml x4	8ml x4	2ml x1	2-8	2year	\$32.00
	Direct Bilirubin*	DBILI	Diazo Salt Method	40ml x4	2ml x1	2ml x1	2-8	2year	\$18.00
	Direct Bilirubin*	DBILI	Vanadate Method	32ml x4	8ml x4	2ml x1	2-8	2year	\$32.00
	Total Protein*	TP	Biuret Method	40ml x4		1ml x1	2-8	2year	\$8.00
	Albumin*	ALB	Bromcresol Green Method	40ml x4		1ml x1	2-8	2year	\$8.00
	Cholinesterase	CHE	Butyrylthiocholine Method	32ml x4	8ml x4		2-8	1year	\$80.00
	Total Bile Acid*	TBA	Enzymatic Cycling Method	32ml x4	8ml x4	2ml x1	2-8	1year	\$350.00
	Renal Function								
	Urea*	UREA	Urease UV Method	32ml x4	8ml x4	2ml x1	2-8	1year	\$24.00
	Creatinine*	CREA	Picric Acid Method	32ml x4	8ml x4	2ml x1	2-30	2year	\$10.00
	Uric Acid*	UA	Uricase-PAP Method	32ml x4	8ml x4	2ml x1	2-8	1year	\$24.00
	Blood Fat								
	Cholesterol*	CHOL	CHOO-PAP Method	32ml x4	8ml x4	2ml x1	2-8	1year	\$24.00
	Triglycerides*	TG	GPO-PAP Method	32ml x4	8ml x4	2ml x1	2-8	1year	\$52.00
	High Density Lipoprotein Cholesterol	HDL-C	Direct Clearance Method	32ml x4	8ml x4	1ml x1	2-8	1year	\$176.00
	Low Density Lipoprotein Cholesterol	LDL-C	Direct Clearance Method	32ml x4	8ml x4	1ml x1	2-8	1year	\$252.00
	Apolipoprotein A 1*	APOA1	Immunoturbidimetric Method	30ml x4	10ml x4	1ml x1	2-8	1year	\$144.00
	Apolipoprotein B*	APOB	Immunoturbidimetric Method	30ml x4	10ml x4	1ml x1	2-8	1year	\$144.00
	Cardiac Muscle								
	Creatine Kinase	CK	IFCC Method	32ml x4	8ml x4		2-8	1year	\$104.00
	Creatine Kinase MB Isoenzyme	CK-MB	AC-Immuno-inhibition Method	32ml x4	8ml x4		2-8	1year	\$336.00
	Lactate Dehydrogenase	LOH	IFCC Method	32ml x4	8ml x4		2-8	1year	\$30.00
	α-Hydroxybutyrate Dehydrogenase	α-HBDH	DGKC Method	32ml x4	8ml x4		2-8	1year	\$30.00
	Pancreas								
	α-Amylase	α-AMY	Gal-G2-CNP Method	40ml x4			2-8	1year	\$432.00
	Leucine Aminopeptidase	LAP	Enzymatic Method	40ml x4			2-8	1year	\$152.00
	Microelements								
	Calcium*	Ca	Arsenazo III Method	40ml x4		2ml x1	2-30	2year	\$12.00
	Chloride*	Cl	Thiocyanate Method	40ml x4		2ml x1	2-30	2year	\$10.00
	Inorganic Phosphorus*	P	Direct UV Method	40ml x4		2ml x1	2-30	2year	\$14.00
	Magnesium*	Mg	Arsenazo III Method	40ml x4		2ml x1	2-30	2year	\$18.00
	Iron*	Fe	Ferrozine Method	32ml x4	8ml x4	2ml x1	2-30	2year	\$58.00
	Saccharides								
	Glucose*	GLU	GOD-PAP Method	32ml x4	8ml x4	2ml x1	2-8	1year	\$16.00
	Glucose*	GLU	Hexokinase Method	32ml x4	8ml x4	2ml x1	2-8	1year	\$28.00
	Fructosamine*	FMN	NBT Method	40ml x4		1ml x1	2-8	1year	\$188.00
	Hemoglobin HBA1C	HBA1C	With Calibrate & Control	One pack					\$890.00

	Specific Protein			4 pcs/pkg		1 pc/pkg	°C	Lifetime	Sales Price For Complete Set
	Immunoglobulin A*	IgA	Immunoturbidimetric Method	30ml x4	10ml x4	1ml x1	2-8	1year	\$212.00
	Immunoglobulin A*	IgA	Immunoturbidimetric Method	15ml x4	5ml x4	1ml x1	2-8	1year	\$106.00
	Immunoglobulin G*	IgG	Immunoturbidimetric Method	30ml x4	10ml x4	1ml x1	2-8	1year	\$212.00
	Immunoglobulin G*	IgG	Immunoturbidimetric Method	15ml x4	5ml x4	1ml x1	2-8	1year	\$106.00
	Immunoglobulin M*	IgM	Immunoturbidimetric Method	30ml x4	10ml x4	1ml x1	2-8	1year	\$212.00

	Immunoglobulin M*	IgM	Immunoturbidimetric Method	15ml x4	5ml x4	1ml x1	2.a	1year	\$106.00
	Complement Component 3*	C3	Immunoturbidimetric Method	30ml x4	10ml x4	1ml x1	2-8	1year	\$154.00
	Complement Component 3*	C3	Immunoturbidimetric Method	15ml x4	5ml x4	1ml x1	2-8	1year	\$78.00
	Complement Component 4*	C4	Immunoturbidimetric Method	30ml x4	10ml x4	1ml x1	2-8	1year	\$154.00
	Complement Component 4*	C4	Immunoturbidimetric Method	15ml x4	5ml x4	1ml x1	2-8	1year	\$78.00
	Prealbumin•	PA	Immunoturbidimetric Method	30ml x4	10ml x4	1ml x1	2-8	1year	\$192.00
	Prealbumin•	PA	Immunoturbidimetric Method	15ml x4	5ml x4	1ml x1	2-8	1year	\$120.00
	Others								
	Lactic Acid*	LA	LA Oxidase Method	32ml x4	8ml x4	2ml x1	2-8	1year	\$78.00
	CRP (C-reactive protein)	CRP		32ml x4	8ml x4		2-8	1year	\$570.00
	HBA1C (Diabetic Test)	HBA1C	142 Test Configuration				2-8	1year	\$680.00
	Calibrator and Control				10 pcs/pkg				
	Multi-Parameter Biochemistry Calibrator	CAL	Calibration		3mlx10		2-8	2year	\$300.00
	Multi-Parameter Biochemistry Control	CON	Contract		3mlx10		2-8	2year	\$220.00
	Detergent				5 pcs/pkg				
	Prothrombin Time	PT	Coagulation Method		100ml x5		2-30	2year	\$76.00
	Highly-effective Detergent (NCH Detergent)	NCH	Detergent		100ml x5		2-30	2year	\$72.00
	Natural Detergent (NCL Detergent)	NCL	Detergent		100ml x5		2-30	2year	\$72.00

COAGULATION REAGENT (SCO-2000) Minimum order 2 sets

No.	Description	Abbreviation	Method	New package					Storage	Price For Complete Set
				R	BUF	C.C'2	OLT	CAL		
				10 pcs/pkg	1 pc/pkg				°C	
	Prothrombin Time	PT	Coagulation Method	2ml x10	30ml x1				2-8	\$95.00
	Activate PartialThromboplastin Time	APTT	Coagulation Method	2ml x10		30ml x1			2-8	\$95.00
	Thrombln Time	TT	Coagulation Method	2ml x10	30ml x1				2-8	95.00
	Fibrinogen	FIB	Coagulation Method	2ml x2	30ml x1		30ml x4	1ml x1	2-8	\$120.00

ELECTROLYTE ANALYZER REAGENT (IO-005) Minimum order 2 sets

No.	Description	Abbreviation	Method	New Package	Storage		Sales Price For Complete Set
					°C	Lifetime	
	Standard Solution A			250ml	2-30	2year	\$60.00
	Standard Solution B			250ml	2-30	2year	\$60.00
	Electrode Reconditioning Solution			20ml	2-30	2year	\$40.00
	Electrode Deproteinization Solution			20ml	2-30	2year	\$40.00
	Electrode Electrolytic Solution			20ml	2 30	2year	\$40.00
	Reference Electrode Solution			20ml	2-30	2year	\$40.00

CLINICAL HEMATOLOGY REAGENT (H-7021)

No.	Description	Abbreviation	Method	New Package	Storage		Sales Price For Complete Set
					°C		
	Diluent (HA-3D)			20L	2-30	2year	\$150.00
	Deprotein (HA-3D)			100ml	2-30	2year	\$116.00
	Hemolysin (HA-3D)			500ml	2-30	2year	\$190.00
	Detergent (HA-3D)			5L	2-30	2year	\$120.00

**One Unit Pack is good for 500 to 600 tests.
Minimum order is \$3,000.00**