



LABOMED, INC.

De-Ionized Water Distiller with Gauge DW-20



Water Processing Volume:	20L/hour
Installation requires:	
a) Water Conductivity:	<400us/cm
b) Water Pressure:	0.1-0.4 Mpa (if no need extra pump)
c) Water Temperature:	5°-45°
Water Speed of Output:	1.5-2L/min
After Purification Conductivity:	>16 Ion < 0.1ppb
Power:	50W
Dimensions:	420mm x 350mm x 500mm
Weight:	25-30 kg