




LABOMED, INC.

www.labomed.com
spectro@labomed.com



LB-122 USB Digital Microscope

Easy to use portable digital microscope with USB output interface.

10× to 200× Magnification.

2.0 Mega Pixels CMOS sensor, interpolate to 5MP with software.

Manual Focus from 10mm to 500mm.

LED illumination with 8pcs LED lamps, brightness adjustable.

Photo management software, Auto White Balance and Exposure.

Powered by USB 2.0 cable, 5V DC.

Metal stand to hold the microscope steady.

Multi-language support: English, Germany, Spanish, French, Russian.

Operating System: Windows Win7, Win 8, Win 10 32 bit & 64 bit,
Mac OS X 10.5 or above.