




LABOMED, INC.

www.labomed.com
spectro@labomed.com



LB-205 Binocular Digital Microscope with Seidentopf Binocular Head, Wide Field Eyepiece, Backward Quadruple Nosepiece, Achromatic Objectives, LED Illumination and 3.2MP Camera and Software

Built-in 3.2 Mega Pixel CMOS camera, vivid images and real color reproduction.

High definition images and videos, high frame rate.

Professional image analyze software supports Windows Vista / Win 7 / Win8 / Mac Operation System.

The software can preview, take photos and video, do image processing and measurement.

Support Multi-language (Arabic, Chinese, English, French, German, Japanese, Polish and etc.).

Compact and flexible, ideally suited for desktop, laboratory worktable.