




LABOMED, INC.

www.labomed.com
spectro@labomed.com



LB-244 Binocular LED Fluorescent Biological Microscope with Wide Field, Epi-Fluorescent Attachment, Seidentopf Binocular Head, Infinite Semi-Plan Achromatic Objectives and Optical System (Infinity Color Corrected System)

Excellent fluorescent image with infinite optical system.

LED illumination with 50,000 hours life span.

No need for cooling system and protection device due to LED cold light source.

465-476nm excited wavelength with high SNR.

Convenient to use, no need warm-up.

More fluorescent filters are available for different applications.