




LABOMED, INC.

www.labomed.com
spectro@labomed.com



LB-353 Digital LCD Stereo Microscope with 8 Inch TFT LCD Screen and 1.3MP Digital Camera

The zoom objective range of optical part is 0.7X~4.5X.

High resolution LCD screen to reverting genuine photo.

With a 1.3MP built-in camera, 8 inches high resolution LCD screen, images can be easily view through the screen.

Can be connected to a PC via USB2.0 interface, images and videos can be stored in the PC.

Cross line, separate screen comparing and freeze functions are available.